



Anti-H₁₁(H₁₁N₂)

Product Name	Anti-H ₁₁ (H ₁₁ N ₂)
Catalog No	IT-003-0015
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	HA ₁ (H ₁₁ N ₂) protein(a.a.19-342) (A/duck/Yangzhou/906/2002 (H11N2)) (Genebank accession # AAY85533)
Applications	Western Blotting (1:200-1:2000)
Specificity	Reacts with HA (H ₁₁ N ₂) protein. Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

