



# Anti-HBsAg (genotype C)

Product Name	Anti- HBsAg (genotype C)
Catalog No	IT-004-001M1
Source	Mouse monoclonal antibody, clone 78C3
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing HBsAg (genotype C) protein (Genebank accession <a href="#">AY040627</a> ).
Applications	ELISA, Neut., IP, IF etc.
Epitope	Conformational
Specificity	Reacts with HBsAg (genotype C). Cross-reactivity to other subtypes not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

