

P2C (HAV)



| | |
|----------------|--|
| Product Name | P2C (HAV) |
| Catalog No | IT-004-051Ep |
| Source | Recombinant protein purified from E.coli |
| Formulation | Each vial contains 100 µg of lyophilized protein in 10mM CBB, pH9.6, 0.1% SDS and 50% glycerol. Reconstitute the protein with 100 µl of Millipore water. |
| Protein | The E.Coli derived recombinant 38.5 kDa protein contains the P2C immunodominant regions, amino acids 1121-1234. |
| Applications | WB standard, antibody ELISA, immunogen, etc. |
| Purity | >90% , based on SDS PAGE |
| Storage | Before reconstitution, stable for 1 year at -20°C from the date of shipment. After reconstitution, stable for a month at 4°C. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |