




p24 (HIV-1/Clade C)

Product Name	p24 (HIV-1/Clade C)
Catalog No	IT-001-019Ep
Source	Recombinant viral protein purified E.coli
Formulation	Each vial contains 100 µg of protein (1 mg/ml) in PBS with 8M urea.
Protein	C-terminal 6xHis tagged HIV-1/clade C p24 protein (Genebank No. AY463234)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	> 95%
Storage	After reconstitution, store at 4°C for short term. For long term storage, aliquot the protein and store at -80°C. Stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: p24 (HIV-1/Clade C)

Note: