

# Anti-M1 (H1N1)



Product Name	Anti-M1 (H1N1)
Catalog No	IT-003-036M2
Source	Mouse monoclonal antibody, clone 36H4.
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG1 in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant protein of M1 from influenza A (A/California/04/2009 (H1N1)) matrix protein (Genbank accession no. ACP44152)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with most M1 from H1N1 viruses. Weak cross-reactivity to other subtypes.
Storage	Store at 4°C; <b>DO NOT FREEZE</b> ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

## ELISA Profile of IT-003-036M2

