




# Anti-GP(Rabies Virus)

---

Product Name	Anti-GP(Rabies Virus)
Catalog No	IT-009-001
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	Ectodomain of GP (Rabies Virus) protein (Genebank No.EU311738)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with GP protein from rabies Virus. Cross-reactivity with glycoprotein from other Rhabdoviruses not tested.
Storage	Store at 4°C; <u>DO NOT FREEZE</u> ; stable for 1 year from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Glycoprotein (Rabies virus; IT-009-001Ep)

Note: