




Gn (Oropouche virus)

Product Name	Gn (Oropouche virus)
Catalog No	IT-017-015Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	6xHis tagged Gn (Oropouche virus) protein (a.a. 24-198)(GenBank No. AJT39490)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: Gn (Oropouche virus)

Note: