

Anti-gp120/gp160(Pan)



Product Name	Anti- gp120/gp160(Pan)
Catalog No	IT-001-001M9
Source	Mouse monoclonal antibody, clone 6B3
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp120(Con_A1/2006) protein (a.a. 34-518)(IT-001-ConA1p)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with HIV-1 Clade A, B, C gp120 proteins as tested. Minimal cross-reactivity to other subtypes.
Storage	Store at 4°C; DO NOT FREEZE ; stable for 6 months from the date of shipment. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-001-001M9

