

# Anti-p24 (HIV-1)



Product Name	Anti-p24(HIV-1)
Catalog No	IT-001-007M14
Source	Mouse monoclonal antibody, clone 80B9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing Gag protein (Clade A) of HIV-1(Genebank No. AY772985)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with p24 and Gag proteins from HIV-1(Clade A)(Clade C).
Storage	Store at 4°C, stable for 6months from the date of shipment. For long term storage, aliquot, and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

