

Anti- 6xHis



| | |
|----------------|--|
| Product Name | Anti- 6xHis tag |
| Catalog No | IT-000-His-Biotin |
| Source | Mouse monoclonal antibody |
| Isotype | IgG1 |
| Formulation | Each vial contains 50 µg biotinylated IgG in 200 µl PBS with 40% glycerol |
| Immunogen | Recombinant protein containing 6xHis tag |
| Applications | ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with 6xHis tag in tagged proteins |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

