

# Anti-Human IgG



---

Product Name	Anti-Human IgG
Catalog No	IT-000-hIg-M2-HRP
Source	Mouse monoclonal antibody, clone 7B2
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing human IgG1 Fc protein
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	HRP labeled monoclonal antibody. Reacts with human IgG and monkey IgG. Minimal cross-reactivity to IgG from others species as tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.