




Anti-gp120/160(HIV-1/Pan)

Product Name	Anti-gp120/160(HIV-1/Pan)
Catalog No	IT-001-0011
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 gp140 proteins (Clade A/B/C)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1(Clade A/B/C) gp120/160 proteins.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
1 : 293 cell expressing gp140(Clade A)
2 : 293 cell expressing gp140(Clade B)
3 : 293 cell expressing gp140(Clade C)

Note: