

Anti-gp120/gp160(Con_A2/2006)



Product Name	Anti-gp120/gp160(Con_A2/2006)
Catalog No	IT-001-001M15
Source	Mouse monoclonal antibody, clone 9G11
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant gp120(Con_A2/2006) protein (a.a. 34-518)(IT-001-ConA2p)
Applications	ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with gp120/gp160 (Con_A2/2006) and gp120 (BL035.W6M.ENV.C1)(A) proteins only as tested. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

ELISA Profile of IT-003-001M15

