




Anti-gp120/160(HIV-1/Clade B)

Product Name	Anti- gp120/160 (Clade B)
Catalog No	IT-001-002
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	HIV-1 gp120 envelope protein ((HXBc2) (Genebank No. K03455)
Applications	Western Blotting (1:200-1:1000)
Specificity	Reacts with HIV-1 gp120/160 (Clade B). Cross-reactivity to gp120/160 from other Clades not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



C: 293 cell extract control
E: 293 cell expressing gp120 (Clade B)

Note: