

# Anti-gp120(89.6P)(HIV-1/Clade B)



---

Product Name	Anti-gp120(89.6P)(HIV-1/Clade B)
Catalog No	IT-001-0022
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium azide.
Immunogen	gp120(89.6P)(HIV-1/Clade B) protein(a.a.34-518) (Genebank accession # <a href="#">U89134</a> ).
Applications	ELISA(1:1000-1:100,000), IP, Neutralizing assay, etc.
Specificity	Reacts with HIV-1 Clade B gp120/160 proteins. Cross-reactivity to other subtype not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

