

# gp120(R2)(HIV-1)

---

Product Name	gp120(R2)(Clade B)
Catalog No	IT-001-0029p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS containing 25% glycerol.
Protein	6xHis tagged HIV-1 gp120(R2)(Clade B) protein (a.a.41-520) (Genebank No. AF128126).
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: gp120(R2) protein

Note: