



hCD317(Δ TMs)

Product Name	hCD317(Δ TMs)
Catalog No	IT-001-317p
Source	Purified protein from transfected 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged ectodomain of human CD317 (Tetherin) protein (a.a.24-162). CD317 has been recently demonstrated to cause retention of retroviral virions on cell surfaces (Nature 451 (7177), 425-430 (2008))
Applications	Interference of Vpr/CD317 interaction, WB standard, , etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: hCD317(Δ TMs)

Note: