



Anti-H1(A/Michigan/45/2015(H1N1))

Product Name	Anti- H1(A/Michigan/45/2015(H1N1))
Catalog No	IT-003-00105M3
Source	Mouse monoclonal antibody, clone 1E9
Isotype	IgG2b
Formulation	Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol
Immunogen	DNA vaccine expressing HA(ΔTM)(H1N1)(A/Michigan/45/2015) protein (IT-003-00105ΔTMp).
Applications	ELISA, IP, IF, WB etc.
Epitope	Linear
Specificity	Reacts with HA(ΔTM)(H1N1)(A/Michigan/45/2015) and some other H1 proteins, minimal or no cross to other HA subtypes tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

