



Anti-H1(A/Michigan/45/2015(H1N1))

| | |
|----------------|--|
| Product Name | Anti- H1(A/Michigan/45/2015(H1N1)) |
| Catalog No | IT-003-00105M4 |
| Source | Mouse monoclonal antibody, clone 3E12 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol |
| Immunogen | DNA vaccine expressing HA(ΔTM)(H1N1)(A/Michigan/45/2015) protein (IT-003-00105ΔTMp). |
| Applications | ELISA, IP, IF etc. |
| Epitope | Conformational |
| Specificity | Reacts with HA(ΔTM)(H1N1)(A/Michigan/45/2015) and some other H1 proteins, minimal or no cross to other HA subtypes tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

