

Anti-H1(H1N1)(A/Michigan/45/2015)



Product Name	Anti-H1(H1N1)(A/Michigan/45/2015)
Catalog No	IT-003-00105m1
Source	Mouse monoclonal antibody, clone 3A10
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1(H1N1)(A/Michigan/45/2015)(IT-003-00105p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reactive to H1(H1N1)(A/Michigan/45/2015) and H1(A/California/07/2009)(H1N1) as tested only. Minimal or no cross to other HA subtypes tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

