

HA(Δ TM)(H1N1) (A/Michigan/45/2015)



Product Name	HA(Δ TM)(H1N1)(A/Michigan/45/2015)
Catalog No	IT-003-00105 Δ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/Michigan/45/2015)(H1N1) (a.a.18-529) (GenBank No. AMV49034).
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq 95\%$
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM)(H1N1)(A/Michigan/45/2015)

Note: