




HA₁ (H1N1)

Product Name	HA ₁ (H1N1)(A/Solomon Islands/3/2006)
Catalog No	IT-003-0011p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA ₁ (H1N1)(A/Solomon Islands/3/2006) protein (a.a.18-344) (GenBank No. ABU99069).
Applications	Western Blotting standard, antibody ELISA, immunogen, etc.
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: purified HA₁ (H1N1) protein