



## HA( $\Delta$ TM)(H<sub>1</sub>N<sub>1</sub>/Solomon Islands/3/2006)

---

Product Name	HA( $\Delta$ TM)(H <sub>1</sub> N <sub>1</sub> /A/Solomon Islands/3/2006)
Catalog No	IT-003-0011 $\Delta$ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA( $\Delta$ TM) (H <sub>1</sub> N <sub>1</sub> /A/Solomon Islands/3/2006)(a.a.18-530)(Genebank No. ABU99069)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

