



HA₁(H₁N₁)(Influenza A)

Product Name	HA ₁ (H ₁ N ₁) (A/Brisbane/59/2007)
Catalog No	IT-003-0012p
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA ₁ (H ₁ N ₁) (A/Brisbane/59/2007) protein (a.a.18-343) (Genebank No. ACA28846)
Applications	WB standard, antibody ELISA, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA₁(H₁N₁) protein

Note: