

# Anti-HA(H1N1(Seasonal))



Product Name	Anti-HA(H1N1(Seasonal))
Catalog No	IT-003-001M13-Biotin
Source	Mouse monoclonal antibody, clone 85F7
Isotype	IgG1
Formulation	Each vial contains 50 µg biotinylated IgG in 100 µl PBS with 40% glycerol
Immunogen	Sequential immunizations using recombinant HA1(H1N1)(A/New Caledonia/20/99(H1N1)) protein (a.a.18-343) (Genebank No. <a href="#">AAP34324</a> ), and HA1(A/California/06/2009(H1N1)) protein (a.a. 18-344)(Genebank accession <a href="#">ACP41935</a> )
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with seasonal H1 including 2009 swine flu HA. Minimal cross-reactivity to other subtypes.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

