



HA(Δ TM)(H₃N₂/Wisconsin/67/X-161/2005)

| | |
|----------------|--|
| Product Name | HA(Δ TM)(H ₃ N ₂ / Wisconsin/67/X-161/2005) |
| Catalog No | IT-003-0041 Δ TMp |
| Source | Viral protein purified from 293 cell culture |
| Formulation | Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS. |
| Protein | C-terminal 6xHis tagged HA(Δ TM) (H ₃ N ₂ / Wisconsin/67/X-161/2005)(a.a.18-530)(Genebank No. ABO37609) |
| Applications | WB standard, antibody ELISA, immunogen, etc |
| Purity | >95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: HA(Δ TM)(H₃N₂/ Wisconsin/67/X-161/2005)

Note: