

HA1(A/Texas/50/2012)(H3N2)



Product Name	HA1(A/Texas/50/2012)(H3N2)
Catalog No	IT-003-00427p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA1(A/Texas/50/2012)(H3N2) (a.a.17-345) (same as HA1 (A/North Carolina/AF2843/2012(H3N2)), Genbank #: AFH57070).
Applications	Protein standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

