




HA(A/Hong Kong/4801/2014)(H3N2)

Product Name	HA(A/Hong Kong/4801/2014)(H3N2)
Catalog No	IT-003-00430ΔTMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 50 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA(A/Hong Kong/4801/2014)(H3N2)(see below for its protein seq)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



Protein SEQ:

qkipgndstatlclghhavpngtiktndrievtnatelvqnssigeicdsphqildgenctlidallgdpqcdgfnkkw
dlfverskaysncypdyvpyaslslvassgtlefnnesfnwtgvtqngtssacirrrssssffsrlnwlthlnykypalnvmt
pnneqfdklyiwgvhpgtdkdqifpyaqssgriivstkrsqavipnigrprirdipsrisiywtivkpgdillinstgnlia
prgyfklrsgkssimrdsapigkcksecitpngsipndkpfqvnritygacpryvkstlkatgmrnvpekqtrgfigaia
gfiengwegmvdgwygfrhqnsgrgqaadlkstqaaidqingklnrligktnekfhqiekefsevegriqdlekyvedtk
idlwsynaellvalenqhtidldsemnklfektkkqlrenaedmngncfkiyhkcdnacigsirngtydhnvyrdealnrr
fqikgvelksygdgsrlvpr