


# HA( $\Delta$ TM)(A/Switzerland/8060/2017(H3N2))



Product Name	HA( $\Delta$ TM)(A/Switzerland/8060/2017(H3N2))
Catalog No	IT-003-00435 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	HA( $\Delta$ TM)(A/Switzerland/8060/2017(H3N2)) (aa 17-529; GISAID# EPI1266285) with 6xHis tagged at C-terminus.
Applications	WB standard, antibody ELISA, Functional analysis, etc
Purity	$\geq$ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)(A/Switzerland/8060/2017(H3N2))

Note: