



HA (B/Malaysia/2506/2004)

| | |
|----------------|--|
| Product Name | HA (B/Malaysia/2506/2004) |
| Catalog No | IT-003-0054p |
| Source | Viral protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | 6xHis tagged HA (B/Malaysia/2506/2004) protein (a.a. 12-547) |
| Applications | Western Blotting standard, antibody ELISA, immunogen, etc. |
| Purity | ≥ 95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: purified HA (B/Malaysia/2506/2004) protein