


# HA<sub>1</sub>(H14N5)(A/mallard/Astrakhan/263/1982)



Product Name	HA <sub>1</sub> (H14N5)(A/mallard/Astrakhan/263/1982)
Catalog No	IT-003-015p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA <sub>1</sub> (H14N5)(A/mallard/Astrakhan/263/1982) protein (a.a.18-347) (Genebank #: <a href="#">P26136</a> )
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA<sub>1</sub>(H14N5)(A/mallard/Astrakhan/263/1982)