



HA(Δ TM) (H₁₄N₅)

Product Name	HA(Δ TM) (H ₁₄ N ₅)
Catalog No	IT-003-015 Δ TMp
Source	Viral protein purified from 293 cell culture
Formulation	Each vial contains 100 μ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged HA(H ₁₄ N ₅) (A/mallard/Astrakhan/263/1982(H14N5)) protein (a.a.18-531) (Genebank No. P26136)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA(Δ TM) (H₁₄N₅) protein

Note: