

# HA( $\Delta$ TM)(B/Massachusetts/2/2012)



---

Product Name	HA( $\Delta$ TM)(B/Massachusetts/2/2012)
Catalog No	IT-003-B19 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged HA(B/Massachusetts/2/2012) (a.a.16-546) (Genbank No. AGL06036)
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

