


# HA( $\Delta$ TM)(B/Wisconsin/01/2010)



---

Product Name	HA( $\Delta$ TM)(B/Wisconsin/01/2010)
Catalog No	IT-003-B7 $\Delta$ TMp
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 $\mu$ g of purified protein (1 mg/ml) in PBS.
Protein	C-terminal 6xHis tagged ectodomain of HA( B/Wisconsin/01/2010) (a.a. 15-548) (Genebank #: EPI271545)
Applications	WB standard, antibody ELISA, immunogen, etc
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA( $\Delta$ TM)(B/Wisconsin/01/2010)