




HA1(A/California/07/2009)(H1N1)

Product Name	HA1(A/California/07/2009)(H1N1)
Catalog No	IT-003-SW12p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS.
Protein	6xHis tagged HA1(A/California/07/2009)(H1N1)(GenBank No. ACP41953; aa 18-345)
Applications	WB standard, antibody ELISA, etc
Purity	≥ 95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. <u>Avoid repeated freezing and thawing</u> . The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: HA1(A/California/07/2009)(H1N1)

Note: