



# Anti-NS5b(HCV)

Product Name	Anti-NS5b (HCV)
Catalog No	IT-004-013M1
Source	Mouse monoclonal antibody
Isotype	IgG1, clone 11B4
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol
Immunogen	Peptide from a conserved region of NS5b polypeptide of Hepatitis C virus (a.a. 2791-2805) (CSSNVSV AHDGAGKR)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with a conserved epitope in NS5b polypeptide of Hepatitis C virus (a.a. 1157-1166) (CSSNVSV AHDGAGKR)
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-003-013M1**

