



gC (Ectodomain)(HSV-2)

| | |
|----------------|--|
| Product Name | gC(Ectodomain)(HSV-2) |
| Catalog No | IT-005-011p |
| Source | Recombinant protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS |
| Protein | C-terminal 6xHis tagged gC(Ectodomain)(HSV-2)(a.a.28-445) (Genebank #: AFM93864). |
| Applications | WB standard, antibody ELISA, immunogen, etc |
| Purity | >95% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

