



# Anti-gG (HSV-2)

Product Name	Anti-gG (HSV-2)
Catalog No	IT-005-013M5
Source	Mouse monoclonal antibody, clone 19D9
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing ectodomain of HSV-2 gG protein (Genebank No. NP_044536) (a.a. 26-310) (IT-005-013p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gG of HSV-2, weakly with HSV-1.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

