

Anti-gH/gL complex (HSV-2)



Product Name	Anti-gH/gL complex (HSV-2)
Catalog No	IT-005-014M12
Source	Mouse monoclonal antibody, clone 38C7
Isotype	IgG2a
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing ectodomain of HSV-2 gH and gL proteins (Genebank No. AFM93846 for gH, and NP_044470 for gL) (a.a. 27-802 and 17-224) (IT-005-014p).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Reacts with gH/gL of HSV-2, but not HSV-1.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

