



# Pp150 (UL32)(Cytomegalovirus)

---

Product Name	Pp150 (UL32)( Cytomegalovirus)
Catalog No	IT-005-021Ep
Source	Viral protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein (1 mg/ml) in 500mM NaCl, 50mM Phosphate buffer (pH 7.4), 200mM imidazole, 8M urea. Reconstitute the protein using 100 µl Ultrapure water.
Protein	C-terminal 6xHis tagged Pp150 (UL32) immunodominant region (a.a. 495-691)
Applications	WB standard, antibody ELISA, etc
Purity	>90%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

