


NS1 (DV4)



| | |
|----------------|--|
| Product Name | NS1 (DV4) |
| Catalog No | IT-006-0055p |
| Source | Viral protein purified from 293 cell culture |
| Formulation | Each vial contains 100 µg of purified protein (1 mg/ml) in PBS. |
| Protein | C-terminal 6xHis tagged Dengue Virus-4 NS1 protein + V5 (Genebank No. ACH61692). |
| Applications | WB standard, antibody ELISA, vaccine, etc |
| Purity | >90% |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |



E: NS1 (DV4) protein

Note: