




NS1 (Zika virus)

Product Name	NS1 (Zika virus)
Catalog No	IT-006-0063Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of protein in PBS with 8M Urea.
Protein	C-terminal 6xHis tagged NS1 (Zika virus, Brazil) antigen (a.a.805-1147) (Genbank Accession No. ALU33341).
Applications	Antibody ELISA, Immunogen, WB standard, etc.
Purity	>90% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: NS1 (Zika virus)