

# Anti-Nucleocapsid protein N(CDV)



Product Name	Anti-Nucleocapsid protein N(CDV)
Catalog No	IT-007-001EM1
Source	Mouse monoclonal antibody, clone 68B5
Isotype	IgG1
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	Recombinant nucleocapsid protein N [Canine distemper virus strain Onderstepoort] protein (a.a. 1-523) (Genbank Accession No. AAK54665) (IT-007-001Ep).
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Recognize nucleocapsid protein N [Canine distemper virus strain Onderstepoort]. Cross-reactivity not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

