

# Anti- Hemagglutinin Protein H (CDV)



Product Name	Anti- Hemagglutinin Protein H (CDV)
Catalog No	IT-007-003M1
Source	Mouse monoclonal antibody, clone 17E10
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing hemagglutinin protein H (Canine distemper virus/strain Onderstepoort)(a.a. 144-604)(Gene Bank Accession No.AAK54669)(IT-007-003p)
Applications	Neut., ELISA, IP, IF etc.
Epitope	Conformational
Specificity	Reacts with H (Canine distemper virus) protein. Cross-reactivity not determined.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

**ELISA Profile of IT-007-003M1**

