

Hemagglutinin Protein H (CDV)



Product Name	Hemagglutinin Protein H (Canine distemper virus)
Catalog No	IT-007-003p
Source	Recombinant protein purified from 293 cell culture
Formulation	Each vial contains 100 µg of purified protein (1 mg/ml) in PBS
Protein	C-terminal 6xHis tagged Hemagglutinin Protein H (Canine distemper virus/strain Onderstepoort) (a.a.144-604)(Gene Bank Accession No. AAK54669)
Applications	Immunogen, antibody ELISA, standard, etc.
Purity	>95%
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80 °C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

