

# M2 Matrix Protein (Rabies virus)



Product Name	M2 Matrix Protein (Rabies virus)
Catalog No	IT-009-005Ep
Source	Recombinant protein purified from E.coli
Formulation	Each vial contains 100 µg of lyophilized protein in PBS with 8M Urea. Reconstitute the protein with 100 µl of Millipore water.
Protein	C-terminal 6xHis tagged M2 Matrix Protein (Rabies virus/strain MRV/genotype 1) (a.a. 1-202)(Genebank Accession No. ABI47939)
Applications	WB standard, antibody ELISA, immunogen, etc.
Purity	>95% , based on SDS PAGE
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.



E: M2 Matrix Protein (Rabies virus)