

Anti-Human RBP-4



Product Name	Anti- Human RBP-4
Catalog No	IT-010-0012
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Peptide from human retinol binding protein 4(hRBP4) (a.a. 178-188) (CLARQYRLIVH).
Applications	WB(1:1000), ELISA(1:2000-1:10000), IP, etc.
Specificity	Reacts with peptide of human RBP-4 (a.a. 178-188) (CLARQYRLIVH).
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.

ELISA Profile of IT-010-0012

