



# Anti- hBNP (preproprotein)

|                |  |
|----------------|--|
| Product Name   | Anti- Human BNP (preproprotein)  |
| Catalog No     | IT-010-002   |
| Source         | Rabbit Polyclonal Antibody   |
| Formulation    | Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.  |
| Immunogen      | Human natriuretic peptide precursor B proprotein protein (a.a.27-134) (Genebank No. <a href="#">NP_002512</a> ).   |
| Applications   | ELISA(1:1000-1:100,000), etc.  |
| Specificity    | Reacts with Human BNP (preproprotein) and mature products  |
| Storage        | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures.   |

**ELISA Profile of IT-010-002**

