

# Anti- Human Apolipoprotein A-I



---

Product Name	Anti- Human APO A1
Catalog No	IT-010-004
Source	Rabbit Polyclonal Antibody
Formulation	Each vial contains 200 µg purified IgG in 100 µl PBS with less than 0.1% gelatin and 0.1% sodium Azide.
Immunogen	Recombinant human APOA1 protein (Genebank accession # X02162). As a major component of the HDL complex (" <b>good cholesterol</b> "), ApoA-I helps to clear cholesterol from arteries.
Applications	ELISA(1:1000-1:100,000), IP, etc.
Specificity	Reacts with human APOA1 protein. Cross-reactivity to other family member not tested.
Storage	Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.