

Anti- Human AMH



| | |
|----------------|--|
| Product Name | Anti- Human AMH |
| Catalog No | IT-010-010M3 |
| Source | Mouse monoclonal antibody, clone 5B6B9 |
| Isotype | IgG1 |
| Formulation | Each vial contains 100µg purified IgG in 100µl PBS with 40% glycerol |
| Immunogen | DNA vaccine expressing Human AMH protein (IT-010-010Ep). |
| Applications | WB, ELISA, IP, IF etc. |
| Epitope | Linear |
| Specificity | Reacts with Human AMH, but not AMH (aa 453-560) as tested. |
| Storage | Keep it at 4°C if used within a month. For long term storage, split it into small aliquots and keep at -80°C. Avoid repeated freezing and thawing. The product will be expired one year after receiving if stored properly. Non-hazardous. No MSDS required. |
| Use Limitation | For research use only, not for use in diagnostic procedures. |

